

Variable Temperature Insert (VTI)

High-Field Air-Core Electromagnet with
Advanced Thermal Management

Controlled Cryogenic Operation.
Stable Vacuum Performance.
Built for Research Applications.

The CRYONANO Cryogenic Dipstick & LN₂ Insert is a modular cryogenic platform designed for precise low-temperature control and electrical measurement in advanced research and industrial environments.

Engineered for operation from liquid nitrogen temperatures (~77 K) to above ambient, the system integrates OFHC copper cold stages, vacuum-compatible electrical feedthroughs, and closed-loop temperature control, enabling stable and repeatable measurements in cryostats, vacuum chambers, and custom experimental setup.

KEY FEATURES

● Thermal Performance

Operates from 77 K to 325 K with rapid thermal equilibration via OFHC copper cold stage, minimized sample temperature gradients, and stable long-duration performance.

● Electrical & Signal Handling

24-pin vacuum feedthrough supports high-voltage lines, sensor signals, thermocouples, and measurement signals while maintaining signal integrity in cryogenic and vacuum environments.

● Mechanical & Vacuum Design

High-vacuum compatible construction with integrated radiation shielding, precision-polished mounting surfaces, and replaceable fixtures for fast sample exchange.

● Temperature Control

Integrated PT1000 sensor, embedded ceramic chip heater, and closed-loop PID regulation support fixed, stepped, and variable temperature modes.



Notes

- Designed for vitrified sample preservation and cryogenic workflows
- Requires proper LN₂ handling and cryogenic safety procedures
- Performance dependent on microscope compatibility and setup

APPLICATIONS

- Low-temperature electrical measurements
- Material characterization at cryogenic temperatures
- Semiconductor device testing
- Sensor calibration and evaluation
- Cryogenic physics experiments
- High-voltage measurements in vacuum

Variable Temperature Insert (VTI)

SPECIFICATIONS

THERMAL & CRYOGENIC PERFORMANCE

Operating Temperature Range	77 K to 325 K (expandable ~70–358 K)
Cooling Method	LN ₂ dip / optional flow channel / cold finger
Cold Stage Material	OFHC Copper
Thermal Conductivity	> 390 W/m·K
Radiation Shielding	Integrated

SAMPLE & MECHANICAL

Sample Size Compatibility	Up to 20 mm × 20 mm (custom supported)
Surface Finish	< 0.1 µm Ra (precision polished)
Mounting	Replaceable fixtures / clamp-based
Flange Options	ISO-KF / CF / custom

TEMPERATURE CONTROL SYSTEM

Sensor Type	PT100 / PT1000 RTD
Sensor Placement	Near sample mount
Heater Type	Embedded ceramic chip heater
Heater Supply	24 V DC ±5%, up to ~2.2 A
Control Method	External PID controller

ELECTRICAL & FEEDTHROUGH

Feedthrough Type	24-pin vacuum feedthrough
Voltage Rating	Up to 500 V DC per pin
Current Rating	~1 A per pin
Signal Types	Power, sensors, thermocouples, test signals
Insulation Resistance	> 10 ⁹ Ω @ 500 VDC
Wiring	Thermally anchored twisted pairs

VACUUM COMPATIBILITY

Sensor Type	PT100 / PT1000 RTD
Sensor Placement	Near sample mount
Heater Type	Embedded ceramic chip heater
Heater Supply	24 V DC ±5%, up to ~2.2 A
Control Method	External PID controller

Cryonano Labs | Head Office: Sector V, WTL Building, BP-5, BP Block, Sector V, Bidhannagar, Kolkata- 700091
Bangalore Office: #999/12, 2nd Floor, GMN Complex, 5th 'A' Cross, 1st Block, HRBR Layout, Babusapalya, Banaswadi,
Outer Ring Road, Bangalore – 560043.
www.cryonano.com | Phone: 097481 81485

Cryonano is not responsible for omissions or errors. Specifications subject to change without notice.
© 2026 CryoNano Instruments Pvt. Ltd. All rights reserved.